

MCL1 Antibody



PACO21565

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

MCL1

Uniprot

Q07820

Synonyms:

TM; EAT; MCL1L; BCL2L3; mcl1/EAT

Applications:

ELISA, WB

Immunogen:

Peptide sequence around aa.148~152(N-N-T-S-T) derived from Human MCL1.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Raji, JK, HL-60 cells using MCL1 Antibody.