

TWIST2 Antibody



PACO21566

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

TWIST2

Uniprot

Q8WVJ9

Synonyms:

DERMO1; MGC117334; bHLHa39


Immunogen:

Peptide sequence around aa.26~30(P-K-R-F-G) derived from Human Twist2.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from A549, Hela and 3T3 cells using Twist2 Antibody.