

ACACA (Ab-79) Antibody



PACO21567

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

ACACA

Uniprot

Q13085

Synonyms:

ACC, ACAC, ACC1, ACCA, ACACAD


Immunogen:

Peptide sequence around aa.78-82(S-M-S-G-L) derived from Human Acetyl-CoA Carboxylase.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa and HUVEC cells using Acetyl-CoA Carboxylase (Ab-79) Antibody.