

CNTFR Antibody



PACO21570

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

CNTFR

Uniprot

P26992

Synonyms:

ciliary neurotrophic factor receptor


Immunogen:

Fusion protein corresponding to a region derived from 23-342 amino acid, of human ciliary neurotrophic factor receptor.

Storage:

Supplied at 0.5mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.3, 0.05% sodium azide and 50% glycerol.

Product Images



Gel: 10%SDS-PAGE Lysate: 40 ug Mouse skeletal muscle tissue lysate
Primary antibody: 1/300 dilution Secondary antibody: Goat anti Rabbit
IgG - H&L (HRP) at 1/10000 dilution Exposure time: 60 seconds.