

AKT1 (Ab-124) Antibody



PACO21580

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:**Protein Background:****Gene ID:**

AKT1

Uniprot

P31749

Synonyms:

AKT; AKT1 kinase; C-AKT; PKB; PKB-alpha

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000


Immunogen:

Synthesized non-phosphopeptide derived from human Akt around the phosphorylation site of serine 124 (S-G-S(p)-P-S).

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from NIH/3T3 cells, treated with PDGF (50mg /ml, 20mins), using Akt (Ab-124) antibody.