

EGFR (Ab-1071) Antibody



PACO21593

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Gene ID:

EGFR

Uniprot

P00533

Synonyms:

EC 2.7.10.1; epidermal growth factor receptor (erythroblastic leukemia viral (v-erb-b) oncogene homolog; avian); Epidermal growth factor receptor precursor; ERBB1


Immunogen:

Synthesized peptide derived from Internal of human EGFR.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from JK cells (Lane 2) and HuvEc cells (Lane 3), using EGFR (Ab-1071) antibody. The lane on the left is treated with synthesized peptide.