

FGFR1/FGFR2 (Ab-463) Antibody



PACO21618

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:**Gene ID:**

FGFR1/FGFR2

Uniprot

P11362/P21802

Synonyms:

Basic fibroblast growth factor receptor 1; bFGF-R-1; BFGFR; CD331; CEK


Immunogen:

Synthesized peptide derived from Internal of human FGFR1/2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from JurKat cells, using FGFR1/2 (Ab-463) antibody.