

BCL6 (Ab-333) Antibody



PACO21621

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:**Gene ID:**

BCL6

Uniprot

P41182

Synonyms:

B-cell CLL/lymphoma 6; B-cell lymphoma 6 protein; B-cell lymphoma 6 protein homolog; Bcl-5; Bcl5


Immunogen:

Synthesized non-phosphopeptide derived from human Bcl-6 around the phosphorylation site of serine 333 (P-Q-S(p)-P-Q).

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 cells, using Bcl-6 (Ab-333) antibody.