

NEUROD1 (Ab-274) Antibody



PACO21628

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:**Gene ID:**

NEUROD1

Uniprot

Q13562

Synonyms:

Neurogenic differentiation factor 1; NEUROD1; NEUROD;


Immunogen:

Synthesized non-phosphopeptide derived from human Neuro D around the phosphorylation site of serine 272 (P-L-S(p)-P-P-).

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa cells, using Neuro D (Ab-274) antibody.