

RPS27A Antibody



PACO21668

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

RPS27A

Uniprot

P62988

Synonyms:

RPS27A; UBA52; UBA80; UBB; UBC

Applications:

ELISA, IHC

Immunogen:

Synthesized peptide derived from Internal of human Ubiquitin.

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using Ubiquitin antibody.

