

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Acts in vitro to stimulate the secretion of adrenocorticotrophic hormone (ACTH). Binds with high affinity to CRF receptor types 1, 2-alpha, and 2-beta. Alison C. Holloway, Am J Physiol Endocrinol Metab, Jul 2002; 283: E165 - E171. Stacie C. Weninger, Endocrinolog.

Gene ID:

UCN

Uniprot

P55089

Synonyms:

UROTENSIN I; UCN;

Immunogen:

Synthesized peptide derived from C-terminal of human urocortin.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using urocortin antibody.

