

ST18 Antibody



PACO21690

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Gene ID:

ST18

Uniprot

O60284

Synonyms:

Suppression of tumorigenicity protein 18; Zinc finger protein 387; ST18; KIAA0535; ZNF387


Immunogen:

Synthesized peptide derived from N-terminal of human ZNF387.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HUVEC cells, using ZNF387 antibody.