
Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Plays a role in nuclear export. Scanlan M. J. , Int. J. Cancer 76:652-658(1998).

Gene ID:

NEMF

Uniprot

O60524

Synonyms:

NY-CO-1; SDCG1; serologically defined colon cancer antigen 1;


Immunogen:

Synthesized peptide derived from internal of human SDCG1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HUVEC cells, using SDCG1 antibody.