

# TADA2A Antibody



PACO21701

---

## Product Information

**Size:**

100ul

**Protein Background:****Reactivity:**

Human, Mouse, Rat

**Gene ID:**

TADA2A

**Source:**

Rabbit

**Uniprot**

O75478

**Isotype:****Synonyms:**

Transcriptional adapter 2-like; ADA2-like protein; TADA2L; KL04P;

**Applications:**

ELISA, WB

**Immunogen:**

Synthesized peptide derived from internal of human ADA2L.

**Recommended dilutions:**


ELISA:1:2000-1:10000, WB:1:500-1:3000

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from COS-7 cells and COLO205 cells, using ADA2L antibody.