

LYL1 Antibody



PACO21729

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:**Protein Background:****Gene ID:**

LYL1

Uniprot

P12980

Synonyms:

Protein lyl-1; Lymphoblastic leukemia-derived sequence 1; LYL1;

Applications:

ELISA, IHC

Immunogen:

Synthesized peptide derived from internal of human Lyl-1.

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using Lyl-1 antibody.

