

HOXB1/HOXA1/HOXD1 Antibody



PACO21732

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:**Protein Background:****Gene ID:**

HOXB1/HOXA1/HOXD1

Uniprot

P14653/P49639/Q9GZZ0

Synonyms:

Homeobox protein Hox-B1; Hox-2l; HOXB1; HOX2l;

Applications:

ELISA, IHC

Immunogen:

Synthesized peptide derived from internal of human HOXA1/B1/D1.


Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Immunohistochemistry analysis of paraffin-embedded human testis tissue, using HOXA1/B1/D1 antibody.