

# HOXB9 Antibody



PACO21738

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

**Protein Background:**

**Gene ID:**

HOXB9

**Uniprot**

P17482

**Synonyms:**

Homeobox protein Hox-B9; Hox-2E; Hox-2.5; HOXB9; HOX2E

**Applications:**

ELISA, WB

**Immunogen:**

Synthesized peptide derived from internal of human HOXB9.

**Recommended dilutions:**


ELISA:1:2000-1:10000, WB:1:500-1:3000

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from K562 cells, using HOXB9 antibody.