

PRKACA Antibody



PACO21739

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

PRKACA

Uniprot

P17612

Synonyms:

EC 2.7.11.11; KAPCA; PKA C-alpha; PKA-alpha; PKACA

Applications:

ELISA, WB

Immunogen:

Synthesized peptide derived from N-terminal of human KAPC A/B.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells and COLO cells, using KAPC A/B antibody.