

EGR1/EGR2 Antibody



PACO21743

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

EGR1/EGR2

Uniprot

P18146/P11161

Synonyms:

AT225; EGR-1; Early growth response protein 1; KROX-24; Krox-24 protein

Applications:

ELISA, WB

Immunogen:

Synthesized peptide derived from internal of human EGR1/2.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa/Jurkat cells, using EGR1/2 antibody.