

APEX1 Antibody



PACO21760

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

APEX1

Uniprot

P27695

Synonyms:

AP endonuclease 1; APE; APE1; APEN; APEX

Applications:

ELISA, WB

Immunogen:

Synthesized peptide derived from internal of human APEX1.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COLO205 cells and HeLa cells, using APEX1 antibody.