

ERCC5 Antibody



PACO21763

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:**Gene ID:**

ERCC5

Uniprot

P28715

Synonyms:

DNA excision repair protein ERCC-5; DNA-repair protein complementing XP-G cells; ERCC5; XPGC; Xeroderma pigmentosum group G complementing protein


Immunogen:

Synthesized peptide derived from internal of human ERCC5.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from colo cells, using ERCC5 antibody.