

CKS2 Antibody



PACO21773

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:**Protein Background:****Gene ID:**

CKS2

Uniprot

P33552

Synonyms:

Cyclin-dependent kinases regulatory subunit 2; CKS-2; CKS2;

Applications:

ELISA, IHC

Immunogen:

Synthesized peptide derived from internal of human CKS2.

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using CKS2 antibody.

