

ATRX Antibody



PACO21782

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:100-1:500

Protein Background:**Gene ID:**

ATRX

Uniprot

P46100

Synonyms:

EC 3.6.1; RAD54L; Transcriptional regulator ATRX; X-linked helicase II; X-linked nuclear protein

Immunogen:

Synthesized peptide derived from N-terminal of human ATRX.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunofluorescence analysis of A549 cells, using ATRX antibody.

