

MAP2K6 Antibody



PACO21800

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:**Gene ID:**

MAP2K6

Uniprot

P52564

Synonyms:

Dual specificity mitogen-activated protein kinase kinase 6; EC 2.7.12.2; MAP kinase kinase 6; MAP2K6; MAPK/ERK kinase 6


Immunogen:

Synthesized peptide derived from N-terminal of human MAP2K6.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 cells, using MAP2K6 antibody.