

LIMK1 Antibody



PACO21805

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:**Protein Background:****Gene ID:**

LIMK1

Uniprot

P53667

Synonyms:

EC 2.7.11.1; KIZ-1; LIM domain kinase 1; LIMK; LIMK-1

Applications:

ELISA, WB

Immunogen:

Synthesized peptide derived from internal of human LIMK1.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from mouse brain cells, using LIMK1 antibody.