

ST5 Antibody



PACO21815

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:**Applications:**

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:**Gene ID:**

ST5

Uniprot

P78524

Synonyms:

Suppression of tumorigenicity 5; HeLa tumor suppression 1; DENN domain-containing protein 2B; ST5; DENND2B


Immunogen:

Synthesized peptide derived from internal of human ST5.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COLO205 cells, using ST5 antibody.