

SP3/SP4 Antibody



PACO21821

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:**Protein Background:****Gene ID:**

SP3/SP4

Uniprot

Q02447/Q02446

Synonyms:

LAG1 longevity assurance; SP260; TMSG1; Tumor metastasis-suppressor gene 1 protein;

Applications:

ELISA, WB

Immunogen:

Synthesized peptide derived from C-terminal of human SP3/4.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells, using SP3/4 antibody.