

NHLH1/NHLH2 Antibody



PACO21822

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

NHLH1/NHLH2

Uniprot

Q02575/Q02577

Synonyms:

Helix-loop-helix protein 1; HEN1; Nescient helix loop helix 1; NSCL-1; NHLH1

Applications:

ELISA, WB

Immunogen:

Synthesized peptide derived from internal of human HEN1/2.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells, using HEN1/2 antibody.