

PFKFB1/PFKFB4 Antibody



PACO21867

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Synthesis and degradation of fructose 2,6-bisphosphate. Sakai A. , J. Biochem. 119:506-511(1996).

Gene ID:

PFKFB1/PFKFB4

Uniprot

P16118/Q16877

Synonyms:

6-phosphofructo-2-kinase/fructose-2; 6-biphosphatase 1; 6PF-2-K/Fru-2; 6-P2ASE liver isozyme; 6-phosphofructo-2-kinase


Immunogen:

Synthesized peptide derived from internal of human PFKFB1/4.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa cells, HepG2 cells, COLO205 cells and 293 cells, using PFKFB1/4 antibody.