

PACO21873

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Gene ID:

PHLDA2

Uniprot

Q53GA4

Synonyms:

Imprinted in placenta and liver protein; Tumor-suppressing subchromosomal transferable fragment candidate gene 3 protein; Tumor-suppressing STF cDNA 3 protein; Beckwith-Wiedemann syndrome chromosomal region 1 candidate gene C protein; p17-Beckwith-Wiedema

Immunogen:

Synthesized peptide derived from internal of human PHLA2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human placenta tissue, using PHLA2 antibody.

