

SFXN4 Antibody



PACO21880

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:**Protein Background:****Gene ID:**

SFXN4

Uniprot

Q6P4A7

Synonyms:

Breast cancer resistance marker 1; Sideroflexin-4; BCRM1; SFXN4;

Applications:

ELISA, WB

Immunogen:

Synthesized peptide derived from N-terminal of human SFXN4.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HUVEC cells, using SFXN4 antibody.