

# CDKL3 Antibody



PACO21892

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Protein Background:**

**Gene ID:**

CDKL3

**Uniprot**

Q8IVW4

**Synonyms:**

Cyclin-dependent kinase-like 3; EC 2.7.11.22; Serine/threonine protein kinase NKIAMRE; CDKL3; NKIAMRE

**Immunogen:**

Synthesized peptide derived from internal of human CDKL3.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunohistochemistry analysis of paraffin-embedded human skeletal muscle tissue, using CDKL3 antibody.

