

PPP1R14C Antibody



PACO21902

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Inhibitor of the PP1 regulatory subunit PPP1CA. Liu QR. , J. Biol. Chem. 277:13312-13320(2002).

Gene ID:

PPP1R14C

Uniprot

Q8TAE6

Synonyms:

Protein phosphatase 1 regulatory subunit 14C; PKCμ; potentiated PP1 inhibitory protein; Kinase-enhanced PP1 inhibitor; Serologically defined breast cancer antigen NY-BR-81; PPP1R14C


Immunogen:

Synthesized peptide derived from internal of human PP14C.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells and COLO205 cells, using PP14C antibody.