

ZNF688 Antibody



PACO21909

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May be involved in transcriptional regulation. Loftus B. J. , Genomics 60:295-308(1999).

Gene ID:

ZNF688

Uniprot

Q8WV14

Synonyms:

Zinc finger protein 688; ZNF688; ZN688;


Immunogen:

Synthesized peptide derived from internal of human ZNF688.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 cells, using ZNF688 antibody.