

DLX4 Antibody



PACO21918

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May play a role in determining the production of hemoglobin S. May act as a repressor. Price J. A. , Hum. Mol. Genet. 7:563-569(1998).

Gene ID:

DLX4

Uniprot

Q92988

Synonyms:

Homeobox protein DLX-4; DLX-7; DLX-8; beta protein 1;


Immunogen:

Synthesized peptide derived from C-terminal of human DLX4.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COLO cells, using DLX4 antibody.