

KAT5 Antibody



PACO21919

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Gene ID:

KAT5

Uniprot

Q92993

Synonyms:

EC 2.3.1.48; EC 2.3.1; 60 kDa Tat interactive protein; Tip60; HIV-1 Tat interactive protein

Immunogen:

Synthesized peptide derived from internal of human TIP60.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using TIP60 antibody.

