

TEP1 Antibody



PACO21937

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Gene ID:

TEP1

Uniprot

Q99973

Synonyms:

Telomerase-associated protein 1; Telomerase protein 1; p240; p80 telomerase homolog;

Immunogen:

Synthesized peptide derived from internal of human TEP1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using TEP1 antibody.

