

IFIH1 Antibody



PACO21945

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Gene ID:

IFIH1

Uniprot

Q9BYX4

Synonyms:

Interferon-induced helicase C domain-containing protein 1; EC 3.6.1; Interferon-induced with helicase C domain protein 1; Helicase with 2 CARD domains; Helicard


Immunogen:

Synthesized peptide derived from C-terminal of human IFIH1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells, using IFIH1 antibody.