

ETNK2 Antibody



PACO21963

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:**Protein Background:****Gene ID:**

ETNK2

Uniprot

Q9NVF9

Synonyms:

Ethanolamine kinase 2; EC 2.7.1.82; EKI 2; Ethanolamine kinase-like protein; ETNK2

Applications:

ELISA, WB

Immunogen:

Synthesized peptide derived from internal of human EKI2.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 cells, using EKI2 antibody.