

# ETV4 Antibody



PACO22004

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Protein Background:**

Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3'.

**Gene ID:**

ETV4

**Uniprot**

P43268

**Synonyms:**

ETS translocation variant 4; Adenovirus E1A enhancer-binding protein; E1A-F; ETV4; E1AF

**Immunogen:**

Synthesized peptide derived from internal of human ETV4.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ETV4 antibody.

