

ACTN1/ACTN2/ACTN3/ACTN4 Antibody



PACO22008

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.

Gene ID:

ACTN1/ACTN2/ACTN3/ACTN4

Uniprot

P12814/P35609/Q08043/O43707

Synonyms:

AAC4; Alpha-actinin 4; F-actin cross linking protein; Non-muscle alpha-actinin 4; actinin


Immunogen:

Synthesized peptide derived from internal of human ACTN4.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 cells, HeLa cells and A549 cells, using ACTN4 antibody.