

# CDK5RAP3 Antibody



PACO22033

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:3000

**Protein Background:**

Potential regulator of CDK5 activity. May be involved in cell proliferation. Regulates CDK5 activity via its interaction with CDK5R1 By similarity.

**Gene ID:**

CDK5RAP3

**Uniprot**

Q96JB5

**Synonyms:**

CDK5 regulatory subunit-associated protein 3; CDK5 activator-binding protein C53; Protein HSF-27; IC53; MSTP016

**Immunogen:**


Synthesized peptide derived from internal of human CDK5RAP3.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from 293 cells, using CDK5RAP3 antibody.