

LAMA1 Antibody



PACO22064

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:100-1:500

Protein Background:

Binding to cells via a high affinity receptor, laminin is thought to mediate the attachment, migration and organization of cells into tissues during embryonic development by interacting with other extracellular matrix components.

Gene ID:

LAMA1

Uniprot

P25391

Synonyms:

LAMA; laminin A; laminin alpha-1;

Immunogen:

Synthesized peptide derived from internal of human LAMA1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunofluorescence analysis of COS-7 cells, using LAMA1 antibody.

