

MP68 Antibody



PACO22096

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Belongs to the small mitochondrial proteolipid family.

Gene ID:

MP68

Uniprot

P56378

Synonyms:

6.8 kDa mitochondrial proteolipid;

Immunogen:

Synthesized peptide derived from internal of human MP68.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma tissue using MP68 antibody.

