

UBA52 Antibody



PACO22106

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Gene ID:

UBA52

Uniprot

P62987

Synonyms:

60S ribosomal protein L40; CEP52;

Immunogen:

Synthesized peptide derived from internal of human RPL40.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using RPL40 antibody.

