

ABHD14B Antibody



PACO22112

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Has hydrolase activity towards p-nitrophenyl butyrate (in vitro). May activate transcription.

Gene ID:

ABHD14B

Uniprot

Q96IU4

Synonyms:

Abhydrolase domain-containing protein 14B; EC 3. -. -. -; CCG1-interacting factor B;

Immunogen:

Synthesized peptide derived from internal of human ABHD14B.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue using ABHD14B antibody.

