
Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Component of the sarcoglycan complex, a subcomplex of the dystrophin-glycoprotein complex which forms a link between the F-actin cytoskeleton and the extracellular matrix.

Gene ID:

SGCA

Uniprot

Q16586

Synonyms:

Alpha-sarcoglycan; Alpha-SG; Adhalin; 50 kDa dystrophin-associated glycoprotein; 50DAG


Immunogen:

Synthesized peptide derived from internal of human SGCA.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells, HepG2 cells and HeLa cells, using SGCA antibody.