

PACO22172

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May be involved in intracellular processing of vitamin B12 (cobalamin). Could play a role in the lysosomal release of vitamin B12 into the cytoplasm.

Gene ID:

ABCD4

Uniprot

O14678

Synonyms:

ATP-binding cassette sub-family D member 4; Peroxisomal membrane protein 69; PMP69; Peroxisomal membrane protein 1-like; PXMP1-L


Immunogen:

Synthesized peptide derived from internal of human ABCD4.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 cells, using ABCD4 antibody.