

DDX51 Antibody



PACO22178

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

ATP-binding RNA helicase involved in the biogenesis of 60S ribosomal subunits By similarity.

Gene ID:

DDX51

Uniprot

Q8N8A6

Synonyms:

ATP-dependent RNA helicase DDX51; DEAD box protein 51; EC 3.6.1. -;


Immunogen:

Synthesized peptide derived from C-terminal of human DDX51.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 cells, using DDX51 antibody.