

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Glycosyltransferase By similarity.

Gene ID:

GCNT6

Uniprot

Q5T4J0

Synonyms:

Beta-1; 3-galactosyl-O-glycosyl-glycoprotein beta-1; 6-N-acetylglucosaminyltransferase 6; EC 2.4.1. -;


Immunogen:

Synthesized peptide derived from internal of human GCNT6.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from K562 cells and NIH/3T3 cells, using GCNT6 antibody.